

# valone

<b>Usage</b>	Upholstery, Bedheads, Bed Valances, Cushions
<b>Width</b>	137cm approx.
<b>Pattern Repeat</b>	n/a
<b>Composition</b>	100% Polyester Back 100% PVC Face
<b>Weight</b>	479grams/m <sup>2</sup> approx.
<b>Flame Retardant</b>	NFPA 260
<b>Abrasion</b>	100,000 Wyzenbeek
<b>Finishes</b>	Resilience - Anti Bacterial, Anti Mold, Bleach Resistant (See reverse side for information on Resilience)
<b>UV Resistance</b>	1000 Hours
<b>Cleaning Care</b>	Mat 9 Prompt cleaning is always recommended. Recommended cleaning steps are listed by type of stain. For best results follow cleaning steps in order.

### Recommended Care and Cleaning Steps

- Step 1:** Mild soap and water, then rinse with fresh water and wipe dry with a clean lint free cloth.
- Step 2:** Wipe affected area with a soft cloth or soft bristle brush using a non-abrasive cleaner. Rinse with fresh water and wipe dry.
- Step 3:** 10% solution of household bleach and water. Rinse with fresh water and wipe dry.
- Step 4:** 25% solution of alcohol and water. Rinse with fresh water and wipe dry.

### Regular Cleaning & Maintenance Step 1 - Dirt, Dust, Grime

Fats & Oils Steps 1 & 2, Food Steps 1 & 2 - Ketchup, Chocolate, Coffee & Tea, Juice

Difficult Stains Steps 1, 2, 3 & 4 - Blood, Urine, Betadine, Eye shadow / Mascara, Lipstick  
Crayon, Grease, Mustard, Ink, Permanent Marker, Paint

<b>Colours Available</b>	Carbon	Nightfall
	Cordovan	Nugget
	Fern	Nutmeg
	Greystone	Pumpkin
	Mudslide	Red Rock

**Lead Time** Stocked in USA, lead time approx. 2 weeks from order unless otherwise quoted.

*Materialised*  
PTY LTD

Sydney 19 Heath Road Blakehurst NSW Australia 2221 p (61) 2 8558 3500 f (61) 2 9546 5402  
sales@materialised.com www.materialised.com  
Adelaide (08) 8271 9185 Brisbane (07) 3393 2188 Hobart (03) 6105 0457 Malaysia (603) 6730 0122  
Melbourne (03) 9815 3033 New Zealand (64) 9 623 0792 Perth (08) 9381 2802 Singapore (65) 6748 9496

specifications

# resilience<sub>sr</sub>

Resilience SR represents the latest technology in Vinyl or Faux Leather protection enhancing performance where it is needed most:-

**Stain Resistance:** with proper care, materials treated with Resilience SR can be easily cleaned and will guard against stains maintaining the desired surface appearance for years to come.

**Anti-Microbial and Anti Bacterial Protection:** Products treated with Resilience SR are engineered to inhibit bacterial growth, as well as fungus, mold and mildew spores, resulting in a more hygienic environment.

**Abrasion Resistance:** Resilience SR provides additional wear properties and surface durability to any Vinyl or Faux Leather, thereby increasing the life of the Material.

Resilience is not visible and does not affect the hand of the materials. Resilience SR is engineered to last, even after repeated cleaning and wear.

Resilience complies with ASTM G 21 / CFFA 120 and AATCC 147 / CFFA 300 inhibiting the growth of the following organisms –

- Aspergillus Niger
- Klebsiella Pneumoniae
- Trichoderma Virens
- Staphylococcus Aureus
- Chaelomium Globosum
- Pseudomonas Aeruginosa
- Penicillium Pinophilum
- Salmonella Choieraesuis
- Aureobasidium Pulluans

Test carried out in USA and available upon request. The information in this cleaning guide refers to performance of resilience in tests conducted under laboratory conditions.

This information does not apply specific guarantees and does not relieve the user from the responsibility of the products safe use and reference cleaning agents.